

réluma^{kiss}

Powered by French Raw Materials & Science



GENÈSE LABS
PARIS



réluma

kiss

A resilient hyaluronic acid filler designed for lip enhancement and perioral rejuvenation

Infused with Vitamin B12, it helps provide a natural rosy hue, enhancing lip color while adding subtle volume, definition, and hydration. Its balanced elasticity and soft texture ensure smooth integration and natural movement, making it ideal for:



Lips
(augmentation & contouring)




Lips

Before




After



Lips

Before



After



Hyaluronic acid concentration	24 mg/mL (Monophasic)
Lidocaine HCl	0.3%
Viscoelasticity	80-160 Pa
Crosslinked agent	BDDE (1,4-Butanediol Diglycidyl Ether)
Injection depth	Subcutaneous
Needle size	30 G * 2

INDIVIDUALIZED TREATMENT

réluma glisse

- Perioral lines (fine lines around the mouth)
- Necklines (horizontal neck wrinkles)
- Décolleté wrinkles
- Crow's feet (fine lines around the eyes)

réluma volumie

- Jawline sculpting
- Chin augmentation
- Cheekbone enhancement
- Deep nasolabial folds

réluma kiss


- Lips (augmentation & contouring)

réluma lush

- Lip augmentation (adding volume and shape)
- Forehead and glabellar fine lines
- Nasolabial folds (smile lines)
- Oral commissures (corners of the mouth)
- Marionette lines

réluma implant

- Jawline contouring (enhances sharpness)
- Chin augmentation (provides structure)
- Chin enhancement (volume and lift)
- Non-surgical rhinoplasty (nose shaping)



GENÈSE LABS
PARIS

WEBSITE www.geneselabs.net | EMAIL info@geneselabs.net
FRANCE HQ 15 Rue des Halles, 75001 Paris, France

